



## XRD-384 COMPATIBLE CARTRIDGES

### Data Sheet

<b>Cat. No.</b>	34-1007	34-1008
<b>Description</b>	8-channel reagent dispenser standard-bore disposable tubing cartridges with plastic-moulded nozzles, compatible with Thermo WellMate Reagent Dispenser	8-channel reagent dispenser small-bore disposable tubing cartridges with plastic-moulded nozzles, compatible with Thermo WellMate Reagent Dispenser
<b>Orifice size</b>	0.023" / 0.58mm	0.015" / 0.38mm
<b>Dispensing accuracy</b>	+/-2% or 1ul	+/-2% or 0.25ul
<b>Dispensing precision</b>	+/-2% or 1ul	+/-2% or 0.25ul
<b>Volume range</b>	20ul and above	1ul and above
<b>Sterile</b>	Yes	Yes
<b>Compatibility</b>	6- to 96-Well Microplates	6- to 96-Well Microplates